Smart System®VR Headset

- Visual Field
 - SMARTracker™ 2 monitors and maintains patient fixation on central target
 - 10-2, 24-2, 30-2 **perimetry**
 - 45 second QuickScreen
 - Full test results in under 3 minutes per eye
 - Supra Threshold
 - Ptosis
 - Esterman (drivers licence)
- Color Vision
- Stereo Testing
- Contrast Sensitivity



glaucoma screening


optimize your patients comfort ...

Clear audio instructions* in the headset prepare your patient optimally for screening. If so preferred, these can also be switched off.

Moreover, Smart System® VR Headset now also offers a comfortable solution for people with reduced mobility.

... while maximizing 'in practice' efficiency!

This testing method will help you determine which patients show an abnormal result and are therefore eligible for further investigations, allowing existing screening equipment to be used more efficiently.



The First Choice in Vision Testing Systems

A HILCO VISION COMPANY

At no time does this product overview serve as medical advice: this authority remains reserved to the ophthalmologist and other authorised medical professionals. The contents of this Simovision® product are based entirely on the specifications available at the time of printing due to the proposed brands and suppliers. Subject to possible printing and typesetting errors.

Simovision® shall not be liable for any damages arising directly or indirectly as a result of consulting this overview or relying on any incorrect information. Products intended for medical professionals only.



more information or a demo?

visit our website

www.simovision.eu

All 'VR Headset' related info gathered on 1 platform:

scientific insights product videos test report & more!





chasing innovation in ophthalmology



premium IOL



retina

cornea

dry eyes



A HILCO VISION COMPANY